

Evac+Chair® 300H Mk4

THE ONLY EVACUATION CHAIR ON THE MARKET WITH A DUAL POSITION SEAT

The Evac+Chair® 300H featuring 400lbs weight capacity and single person operation.



Features

- 1 Adjustable patient head restraint**
Limits movement, securing the passenger in place
- 2 Friction tracks (Evac+Track technology)**
The self braking reinforced tracks control speed of descent for easy single person operation
- 3 Ergonomic adjustable handle**
For ease of deployment and compact storage
- 4 Reflective decals**
For easy location in situations of poor visibility
- 5 Quick release safety belts (Head, Chest, Thigh)**
For increased passenger safety
- 6 Heavy duty wheels**
Ideal for operation over uneven terrain
- 7 Lockable rear wheels**
Can be applied in situations where the chair needs to remain in a stationary position
- 8 Dual position seat**
Dual position seat allows the user to decide each time which seating position they desire



Includes

- LIFETIME WARRANTY
- 3 Passenger Safety Belts (Head, Chest, Thigh)
- Photoluminescent Signage
- Wall Mounting Brackets
- Dust Cover
- User Guide
- Instructional DVD

Benefits

- Designed in the US
- FDA Registered
- Provides a quick, easy and safe emergency evacuation
- Easy to use
- One person operation
- Lightweight
- No heavy lifting required
- Enables small attendants to easily move passengers three times their size
- Easy to open, ready to use in seconds
- Folds for compact storage
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act
- Wall mountable



Flat seat

with zipped side panels and over thigh quick release passenger restraint



Hammock seat

simply release zipped side panels, fasten together on underside, retaining use of over-thigh quick release passenger restraint



Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	41"	20"	8"	19lbs	400lbs

VIDEO LINK: evac-chair.ca/watch-video/

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®. Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11, EXCLUSIVE LIFETIME WARRANTY.

For further information please contact:

Rescue 7 Inc.

T 905-474-0770 E sales@rescue7.net

W www.rescue7.net

EVAC+CHAIR®
CANADA

