

Evac+Chair® 400H

THE WORLD'S BEST SELLING EVACUATION CHAIR JUST GOT EVEN BETTER!

The Evac+Chair® 400H featuring a weight capacity of 440lbs and single person operation.

MODEL BENEFITS

- Reinforced Adjustable Friction Tracks
- 440lbs weight capacity
- 5 Point Safety Harness
- Passenger Armrests and Footrest



FEATURES

- 1. Adjustable patient head restraint**
Limits movement, securing the passenger in place
- 2. Adjustable friction tracks**
(EVAC+TRAC technology)
Adjustable self braking reinforced tracks allow for tension change to control speed of descent to suit most staircases
- 3. Ergonomic adjustable handle**
For ease of deployment and compact storage
- 4. Passenger armrests**
For added passenger comfort and safety.
- 5. 5 Point quick release safety harness**
For increased passenger safety
- 6. Heavy duty wheels**
Ideal for operation over uneven terrain
- 7. Large lockable rear wheels**
Can be applied in situations where the chair needs to remain in a stationary position
- 8. Passenger footrest**
Allows for stable passenger seating position during stairway operation
- 9. Dual position adjustable seat**
Adjustable seat allows the passenger to choose desired seating position

Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	41"	20"	8"	33lbs	440lbs

BENEFITS

- Designed in the US
- FDA Registered
- Provides a quick, easy and safe emergency evacuation
- Easy to use
- One person operation
- Lightweight
- No heavy lifting required
- Enables small attendants to easily move passengers three times their size
- Easy to open, ready to use in seconds
- Folds for compact storage
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act
- Wall mountable

INCLUDES

- LIFETIME WARRANTY**
- 5 Point Quick Release Harness
- Photoluminescent Signage
- Wall Mounting Brackets
- Dust Cover
- User Guide

VIDEO LINK:

evac-chair.com/watch-video/

For further information please contact:

EVAC+CHAIR®
CANADA

