

THE WORLD'S NO.1 STAIRWAY EVACUATION CHAIR

Evac+Chair® Power 800


THIS ELECTRIC/BATTERY POWERED STAIR CLIMBER ALLOWS FOR EASY AND SAFE STAIR ASCENT & DESCENT

Can be used in emergency evacuations, power outages, and daily patient transport to safely and easily move a passenger up and down a flight of stairs.



Features

- 1 Height adjustable head cushion**
Headrest with padded cushion for maximum comfort for the passenger during transit
- 2 Easy to use control panel**
3 Speed control panel with on/off power and ascent and descent buttons
- 3 Ergonomic adjustable handle**
For ease of deployment and compact storage
- 4 Additional operator handles**
For added stability when using additional operators during transportation
- 5 Extra long quick release four point harness**
For increased passenger safety
- 6 Rechargeable lithium Ion battery pack**
Fully charged in 1.5 hours
(Equivalent to Going Up and Down a 5 Story building flight of stairs 30 times)
- 7 High strength aluminium frame**
Increased strength and robustness for continued use
- 8 LED safety light and photoluminescent branding**
For ease of use in situations of poor visibility
- 9 Adjustable stair track**
Tension control system for perfect transportation with tool included
- 10 Lockable rear wheels**
Can be applied in situations where the chair needs to remain in a stationary position



Includes

- Battery Charger (US/EU/UK)
- Rechargeable Battery
- 2 Keys
- User Guide
- Passenger Headrest
- Passenger Safety Straps
- Photoluminescent Wall Sign
- Wall Hooks
- Dust Cover
- Manufacturers Warranty

Benefits

- FDA Registered
- CE & ISO 9001 Certificated
- Robust and compact in size
- Provides a quick, easy and safe emergency evacuation
- Easy to use
- One person operation
- No heavy lifting required
- Enables small attendants to easily move passengers three times their size
- Easy to open, ready to use in seconds
- Folds for compact storage
- Flame retardant materials
- Wipeable materials and anti-microbial paint for infection control
- Exceeds fire safety requirements and emergency response guidelines established to support the Americans with Disabilities Act



Dimensions	Height	Width	Depth	Weight	Weight Capacity
Folded	45"	19.5"	10.6"	66lbs	350lbs

VIDEO LINK: evac-chair.ca/watch-video/

For over 35 years EVAC+CHAIR® is America's Original and World's Number One Selling Evacuation Chair. Don't Risk Lives, Beware of Poor Imitations and Inferior Substitutions. There is only one EVAC+CHAIR®. Unsurpassed Quality and Workmanship, Successful in Saving Lives in 9/11.

For further information please contact:

Rescue 7 Inc.

T 905-474-0770 E sales@rescue7.net

W www.rescue7.net

EVAC+CHAIR®
CANADA

